

CCMC-T05

5 metre camera cable for PMW-10MD



Overview

Extends flexibility of medical camera

The CCMC-T05 camera cable allows you to extend the PMW-10MD medical two-piece camera. The PMW-10MD has a compact camera head unit, that can be attached easily either to a surgical microscope, a surgical boom arm, or to a ceiling, and a separated camera control unit. The 5 m cable increases the distance between the camera head and the camera control unit.

- **MDD approved**

The camera cable meets MDD approval in the EU.

Features

- **Extends distance between camera head and camera control unit**

The 5 m camera cable is fully MMD approved in the EU and extends the flexibility of the PMW-10MD in surgical conditions.

Compatible Products

Medical Imaging Cameras



PMW-10MD

Full HD Medical 2-piece
Camera